



## 853 Antares

Cervical collar set 3 pieces with carrying case  
Art. 844 - 846 - 848



### FEATURES

ANTARES cervical collar is an highly professional instrument indispensable for the emergency operator. It is composed by 2 pieces to be put on the front and back part of the neck enabling thus the perfect immobilization of the rachis. Central hole for tracheotomy and easy positioning thanks to the rip-off fastening which different colours indicate the various sizes. It is made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different sizes.  
Art. 844 - 846 - 848

### MAIN FEATURES

|        |                                                                                     |
|--------|-------------------------------------------------------------------------------------|
| Usage  | Immobilization                                                                      |
| Trauma |  |

---


SEE ALSO

---



**Sirius 874**

Cervical collar set 4 pieces with carrying case Art. 862 - 864 - 866 - 868



**Sirius 876**

Cervical collar set 6 pieces with carrying case



**Antares 844**

Rigid cervical collar in 2 pieces size STOUT



**Antares 846**

Rigid cervical collar in 2 pieces size STANDARD



**Antares 848**

Rigid cervical collar in 2 pieces size MEDIUM

---