



11064

Manual suction unit



FEATURES

HAND SUCTION UNIT

Designed to provide quick and effective pharyngeal and tracheal suction for emergency use.

The adjustable stroke knob provides different suction pressure.

Supplied with 250 ml bottle and 2 non sterile adult cannulas 17 cm (Ø 6 mm) and 19 cm (Ø 10 mm).

Max vacuum: -400 mmHg
Peak free airflow: >20 l/min
External diameter of catheter connection: 10 and 6 mm
Weight: 260 g
Size: 230x60x190 mm

MAIN FEATURES

Length (cm)	23
Height (cm)	19
Width (cm)	6
Weight (kg)	0,26

ACCESSORY OF



1366

Red knapsack "115" with side
pockets and compartments
COMPLETE
