



## 1362

“Ambumed” rescue knapsack



### FEATURES

New modular knapsack, practical, compact, simple and suitable for rescue volunteers and professionals. Externally it is made of blue tear-resistant double sewing nylon, sturdy 10 mm zips, reflective high visibility stripes on all sides and a closed- cell PVC bottom, with plastic supports. On the inner side there is a protective hard shell structure made of thermoplastic material, which consists of a central wide space, characterized by a double-zip opening, and 4 external side pockets. The central space is supplied with Velcro® stripes for holding 4 coloured pockets, each ones supplied with a large transparent window reinforced with double sewing, and 2 professional bags to store and organize different tools and materials. This knapsack is completed by a frontal pocket with inner straps to secure different objects, 3 anti-slip rubber handles (1 on the upper side and 2 on lateral sides) to carry by hand and, on the back side, 2 shoulder bands with reflecting stripes.

Ambumed is complete with 4 pockets (Art. 1362-002) and 2 professional bags (art. 1362-001).

**Dimensions: 53 x 53 x 25 cm (LxHxD).**

### MAIN FEATURES

|             |       |
|-------------|-------|
| Width (cm)  | 53    |
| Height (cm) | 53    |
| Depth (cm.) | 25    |
| Weight (kg) | 4,1   |
| Material    | Nylon |

IMAGES



SEE ALSO



**1325/RB**  
Compact professional rescue  
backpack Antartide



**Rescuemed 1343**  
Rescue knapsack



**1354**  
Red knapsack "115" with side  
pockets and compartments

ACCESSORIES



**1362-001/R**  
Red squared inside bag for  
knapsacks



**1362-001/G**  
Yellow squared inside bag for  
knapsacks



**1362-001/V**  
Green squared inside bag for  
knapsacks



**1362-002/A**  
Orange rectangular inside bag for  
knapsacks



**1362-002/B**  
Blue rectangular inside bag for  
knapsacks



**1362-002/R**  
Red rectangular inside bag for  
knapsacks