



7240/INC PROOF Frog

"FROG PLUS INC PROOF" self loading trolley for incubators
predisposed for certified fixing system



FEATURES

Made with predisposition for our certified fixing device Art. 955 Proof, it is perfect for transport of incubators on rescue vehicles. Designed without back rest, side protections and rigid platform. This stretcher is part of "FROG PLUS" self loading stretchers range.

Device certified according to the European harmonized safety standards UNI EN1865 - UNI EN 1789.

MAIN FEATURES

Pediatric	
Height (cm)	82
Max loading height (cm)	66
Min loading height (cm)	60
Folded height (cm)	27
Width (cm)	56
Length (cm)	195
Length with folded trolley (cm)	166
Weight (kg)	38,4
Loading capacity (kg)	250

SEE ALSO



Frog 7240 PROOF
"FROG PLUS" self loading stretcher
with trendelenburg certified

ACCESSORIES



12050
Extended trolley incubator
interface



12048
Certified interface for incubator
fixation for self loading trolleys



955 PROOF
EN 1789 fixing system for Proof roll-
in stretchers
