



8506

Pulse oximeter Oxy-110



FEATURES

- Large 3.5 inch auto rotation TFT display for SpO2, PI, Plethysmogram, PR and temperature
- patented MoveOxy™ SpO2 technology with artifact removal and anti-motion performance
- various choice of SpO2 probes for adult, pediatric and neonate
- built-in lithium battery up to 18 hours continuous monitoring
- adjustable alarm limit, pitch tone level, key tone and brightness
- up to 580 hours SpO2 data storage
- upgradeable to infrared ear temperature

Supplied with:

- cover
- adult clip probe
- charging base
- USB cable

MAIN FEATURES

Thickness (cm)	2,5
Height (cm)	15,8
Width (cm)	7,3
Weight (kg)	0,23
Cardio	

IMAGES



SEE ALSO



11220 Oxytest
OXYTEST portable pulse oximeter
complete with adult sensor, USB
cable and AAA battery



11200 Vital Test
Handheld pulse oxymeter
complete with adult sensor

ACCESSORIES



8506-001
Sensore SpO2 adulto riutilizzabile
per pulsiossimetro OXY-110



8506-002
Sensore SpO2 pediatrico
riutilizzabile per pulsiossimetro
Oxy-110



8506-003
Sensore SpO2 neonatale
riutilizzabile per pulsiossimetro
Oxy-110
