




19146-956

"MX-2-956" double fixed space stretcher support with ramp of mm. 350, equipped with certified fixing system art. 956 - EN1789 certified



FEATURES

Double space fixed stretcher support with ramp mm. 350. On request supplied with UNI EN 1789 certified blocking devices Art. 956.
Gas spring included.

MAIN FEATURES

Width (cm)	63
Length (cm)	223
Weight (kg)	68,2
Loading capacity (kg)	400
Height (cm)	21,5

ACCESSORIES



19110

Stainless steel plate mm 1940 x 572
for stretcher supports



787L

Gas spring for ramp

STANDARD ACCESSORIES



821G

Gas spring for ramp
