



## 658/Y Quick

Folding sedan chair with 2 wheels – yellow/black



### FEATURES

Designed to safely carry patient, and to be opened/closed in simple and safe way. A special mechanism avoids accidental closure of the chair. Supplied with 2 wheels Ø 125 mm, 2 rear levers and 2 front extractable ones. Easy to use. Back rest and seat made of thermoplastic material easy to clean. Supplied complete with black bag containing 3 belts (Art. 590/N).

Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789.

### MAIN FEATURES

10G Certification	
Foldable	
Wheel type	2 Small Fixed Wheels
Type	
Height (cm)	90
Width (cm)	51
Length (cm)	71
Max length (cm)	125
Epaisseur pliée (cm)	23
Weight (kg)	10
Loading capacity (kg)	170
Folded height (cm)	90
Cardio	
Material	Painted steel

---

## SEE ALSO

---



**Quick 659/Y**

Folding sedan chair with 4 wheels,  
footboards and arms - yellow/black

---

---

## ACCESSORIES

---



**618**

Fixing system for sedan chairs



**612/AN-2/MEB**

Adjustable thoracic abdominal belt  
orange/black with MeBer band and  
4 points fastening for self loading  
stretchers and sedan chairs



**612/N-2/MEB**

Adjustable thoracic abdominal belt  
black with MeBer band with 4  
points fastening for self loading  
stretchers and sedan chairs



**590/N**

Black bag with 3 belts for sedan  
chairs

---

---

## STANDARD ACCESSORIES

---



**590/N**

Black bag with 3 belts for sedan  
chairs

---