



## 659/Y Quick

Folding sedan chair with 4 wheels, footboards and arms - yellow/black



### FEATURES

Indispensable device for rescue operator. Supplied with 4 wheels (2 rear fixed Ø 125 mm and 2 front rotating with brake Ø 100 mm), foldable arms and footrest, 4 rear handles and 2 front extractables ones. This chair grants safety and comfort to the patient. Back rest and seat made with thermoplastic material easy to clean. Supplied complete with black bag containing 3 belts (Art. 590/N).

Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

### MAIN FEATURES

10G Certification	
Foldable	
Footboards and arms	
Wheel type	2 Rotating Wheels + 2 Small Fixed Wheels
Type	
Height (cm)	90
Width (cm)	51
Length (cm)	71
Max length (cm)	126
Epaisseur pliée (cm)	23
Weight (kg)	14
Loading capacity (kg)	170
Folded height (cm)	99
Cardio	
Seat length (cm)	47
Material	Painted steel

---

## SEE ALSO

---



**Quick 658/Y**  
Folding sedan chair with 2 wheels –  
yellow/black

---

---

## ACCESSORIES

---



**618**  
Fixing system for sedan chairs



**612/AN-2/MEB**  
Adjustable thoracic abdominal belt  
orange/black with MeBer band and  
4 points fastening for self loading  
stretchers and sedan chairs



**612/N-2/MEB**  
Adjustable thoracic abdominal belt  
black with MeBer band with 4  
points fastening for self loading  
stretchers and sedan chairs



**590/N**  
Black bag with 3 belts for sedan  
chairs

---

---

## STANDARD ACCESSORIES

---



**590/N**  
Black bag with 3 belts for sedan  
chairs

---