



## 939

Instrument support for stretchers



### FEATURES

Instrument support to be fixed on Me.Ber. stretcher frames' to transport medical equipments, defibrillators and monitors. It is supplied with belts to fix devices and is foldable to be stored inside ambulances.

Dimensions: open cm 34x56x37h, closed height cm 8

### MAIN FEATURES

Width (cm)	56
Depth (cm.)	34
Height (cm)	37
Epaisseur pliée (cm)	8
Weight (kg)	5,6
Loading capacity (kg)	30

### IMAGES



---

## SEE ALSO

---



**Mercury Lite 7082/4RG Proof**  
Certified aluminium variable height  
stretcher high version with 4  
rotating wheels

---

---

## ACCESSORY OF

---



**Mercury Cinque 7075/4RG PROOF**  
Certified 5 levels self-loading  
stretcher with 4 swivel wheels



**Mercury Cinque 7077/4RG PROOF**  
Certified 5 levels self-loading  
stretcher with 4 swivel wheels -  
high version



**Mercury Lite Cinque 7095/4RG  
Proof**  
Certified aluminium 5 levels self-  
loading stretcher with 4 swivel  
wheels



**Mercury Lite Cinque 7097/4RG  
Proof**  
Certified aluminium 5 levels self-  
loading stretcher with 4 swivel  
wheels - high version

---



**Mercury EL 7092 PROOF**  
Certified self-loading stretcher



**Mercury EL 7094 PROOF**  
Certified self-loading stretcher -  
high version



**7265 Proof**  
Aluminium self loading stretcher  
"Frog lite Advance" with  
extractable handles and  
Trendelenburg

---